

## Instrument Decontamination Bath 3L (includes basket, drip tray and lid)

**Product code: MJ-DSS2711**

Dimensions:  
257x180x128mm

Volume: 3 litre

Material: Polypropylene

Colour: Blue



Used for soaking and cleaning of instruments and devices.

- Lightweight, strong and durable – 50% lighter than stainless steel.
- Designed to withstand high impact – will not dent, rust or strain.

Instructions: To decontaminate - cover instruments with **Empower Foam** or **Enzymax Spray Gel** or soak and wash in a **ProDet** clinical detergent solution. The cleaning/soaking container can be autoclaved.



Empower Foam 24oz bottle ea

KE-10-4224S



Enzymax Spray Gel 709ml

HF-IMS1229



ProDet Clinical Detergent 5 Litre  
Container

MJ-PD5