

For lifelike dentures with multi-layered shape design.



> FACTS & FEATURES

What? • lifelike, detailed premium anterior teeth made of MRP composite

- Benefits?**
- **ideal integration in the remaining natural dentition**
through detailed, multi-layered moulds
 - **lively play of light**
through elaborate surface texture
 - **good shade integration with natural remaining dentition**
thanks to high shade accuracy to the shade guide
 - **reliable masking of secondary/tertiary structures**
thanks to high (basal) tooth volume
 - **excellent durability**
thanks to high load capacity and abrasion stability

What for? • ideal solution for partial and hybrid dentures



> BASIC INFORMATION

Characteristics

Mould

- natural three-dimensional tooth geometry for optimum integration in the remaining dentition
- lifelike, distinctive incisal characteristics and natural labial curvature
- broad tooth neck for ideal interdental space/papilla modeling
- high (basal) volume for masking secondary/tertiary structures

Special features

- excellent chromaticity for good masking of tertiary structures
- elaborate surface structure for lively play of light
- 3D-layered structure for natural shade gradient

Recommended indications

Ideal prosthetic solution for:

- Partial prostheses
- model casting, attachment and bar constructions
- Telescope and cone constructions

Occlusal/setup concept

★ Especially suited for:

BLP (BIO-Logical Prosthetic) according to Dr. E. End

Available shades

• 28 x VITA SYSTEM 3D-MASTER® shades

0M1, 0M3, 1M1–5M3

• 15 x VITA classical A1–D4® shades

A1–D4 (without B1)

Available moulds/sets

- 21 x upper anterior tooth moulds (divided into four shape groups)
- 8 x lower anterior tooth moulds

> OVERVIEW OF SHAPES

Upper Anteriors							
 ovoid	O1S*  10.8 7.6	O2L  11.5 8.4	O3M  10.8 8.6	O4L  12.7 9.0	O5L  11.4 9.1	O6L  13.1 9.7	
	T6S**  10.4 7.4	T2S*  10.1 8.3	T1S  10.2 8.7	T3M  11.7 8.4	T4M*  11.3 8.8	T5M  11.1 8.8	T6L  12.4 8.9
 triangular	T7M  11.6 8.8	T8L*  11.6 9.2	T9L  13.1 9.4				
	X1M*  11.0 8.8	X2L  11.6 8.7	X3L  11.8 9.0				
 rectangular	Z1S  9.8 8.1	Z2S  10.1 8.7					
 square							
Lower Anteriors							
	L1S*  8.8 4.9	L2M*  10.4 5.2	L3M*  10.5 5.2	L4M*  9.4 5.4	L5M*  10.4 5.5	L6L  11.4 5.8	L7L  10.4 5.5
	L8L  13.0 6.3						

*) Bleached shade(s) OM1 and OM3 exclusively for selected mould(s).

**) Mould(s) only available in VITA SYSTEM 3D-MASTER shade(s).

VITA PHYSIODENS®

RECOMMENDED COMBINATIONS





Combinations according to indications:

For partial and hybrid dentures

- VITA PHYSIODENS Anterior with VITA PHYSIODENS Posterior

VITA PHYSIODENS Anterior		VITA PHYSIODENS Posterior	
UJ / No.  mm	LJ / No.  mm	UJ / No.  mm	LJ / No.  mm
01S 42.3	L1S 33.0	21E 30.5	21E 32.7
02L 45.7	L2M 34.5	22E 32.9	22E 33.7
03M 47.2	L3M 36.3	23E 33.5	23E 34.6
04L 49.2	L6L 38.2	24E 34.8	24E 36.1
05L 50.1	L6L 38.2	24E 34.8	24E 36.1
06L 51.3	L8L 43.5	25E 36.7	25E 37.3

T6S 41.0	L1S 33.0	21E 30.5	21E 32.7
T2S 44.9	L2M 34.5	21E 30.5	21E 32.7
T1S 45.3	L1S 33.0	20E 29.4	20E 31.0
T3M 45.6	L3M 36.3	22E 32.9	22E 33.7
T4M 45.6	L4M 35.5	22E 32.9	22E 33.7

VITA PHYSIODENS Anterior		VITA PHYSIODENS Posterior	
UJ / No.  mm	LJ / No.  mm	UJ / No.  mm	LJ / No.  mm
T5M 46.2	L3M 36.3	23E 33.5	23E 34.6
T6L 47.3	L5M 37.0	23E 33.5	23E 34.6
T7M 47.5	L7L 38.2	23E 33.5	23E 34.6
T8L 49.0	L5M 37.0	24E 34.8	24E 36.1
T9L 49.6	L6L 38.2	25E 36.7	25E 37.3

X1M 46.3	L3M 36.3	23E 33.5	23E 34.6
X2L 46.2	L4M 35.5	23E 33.5	23E 34.6
X3L 47.8	L5M 37.0	24E 34.8	24E 36.1

Z1S 44.1	L1S 33.0	21E 30.5	21E 32.7
Z2S 47.5	L5M 37.0	22E 32.9	22E 33.7

upper jaw (UJ) / lower jaw (LJ)

The information contained in the combination table is only intended as a guide and can vary, depending on the patient's individual situation.

› ASSORTMENTS AND STORAGE



VITA PHYSIODENS mould chart

Handy metal case with all moulds for tooth selection:

- VITA PHYSIODENS Anterior and Posterior



VITA PHYSIODENS tooth storage

Anterior:


- X-small*: 67 Sets / 402 Anterior
- Small*: 134 Sets / 804 Anterior
- Medium*: 268 Sets / 1,608 Anterior
- Large*: 400 Sets / 2,400 Anterior

Available in VITA classical A1–D4® and VITA SYSTEM 3D-MASTER® shades.**

*) Note: You will receive the most frequently used moulds and shades with the packaged tooth storage. We are also happy to put together storage systems according to your individual specifications.

**) see available shades

› SYSTEM COMPONENTS





**PREMIUM
PROSTHETIC**



reliable masking

three-dimensional,
natural moulds